

Human KLRB1/CD161 protein (Recombinant) (C-hFc) (STJP019893)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3820
Gene Symbol	KLRB1
UniProt ID	KLRB1_HUMAN
Immunogen Region	Gln67-Ser225
Immunogen Sequence	QKSSIEKCSVDIQQSRNKTT ERPGLLNCPYWQQLREKCL LFSHTVNPWNNSLADCSTKE SSLLLIIRDKDELIHTQNLIR DKAILFWIGLNFSLSEKNWK WINGSFLNSNDLEIRGDAKE NSCISISQTSVYSEYCSTEI RWICQKELTPVRNKVYPDS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---