

Human c-Kit/CD117 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019897)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	3815
Gene Symbol	KIT
UniProt ID	KIT_HUMAN
Immunogen Region	Gln26-Pro520
Immunogen Sequence	QPSVSPGEPSPPSIHPGKSD LIVRVGDEIRLLCTDPGFVK WTFEILDETNNKQNEWITE KAEATNTGKYTCTNKHGLSN SIYVFVRDPAKLFLVDRSLY GKEDNDTLVRCPLTDPEVTN YSLKGCQGKPLPKDLRFIPD PKAGIMIKSVKRAYHRLCLH CSVDQEGKSVLSEKFILKVR PAFKAVPVVSVSKASYLLRE GEEFTVTCTIKDVSSSVYST WKRENSQTKLQEKYNSWHH
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.