

Human IL-15RA/CD215 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019899)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3601
Gene Symbol	IL15RA
UniProt ID	I15RA_HUMAN
Immunogen Region	Ile31-Thr205
Immunogen Sequence	ITCPPPMSVEHADIWVKSYS LYSRERYICNSGFKRKAGTS SLTECVLNKATNVAHWTTPS LKCIRDPALVHQRPAAPPSTV TTAGVTPQPESLSPSGKEPA ASSPSSNNTAATTAAIVPGS QLMPSPKSPSTGTTEISSHES SHGTPSQTTAKNWELTASAS HQPPGVYPQGHSDDT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---