

Human IL-12 R beta 1/CD212 protein (Recombinant) (C-hFc) (STJP019901)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as 0.22 Mu m filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -70 °C. This product is stable at ≤-70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3594
Gene Symbol	IL12RB1
UniProt ID	I12R1_HUMAN
Immunogen Region	Cys24-Glu540
Immunogen Sequence	CRTSECCFQPPYPDADSGS ASGPRDLRCYRISSDRYECS WQYEGPTAGVSHFLRCCLSS GRCCYFAAGSATRLQFSDQA GVSPLYTVTLWVESWARNQT EKSPEVTLQLYNSVKYEPL GDKVSKLAGQLRMWETPD NQVGAEVQFRHRTSPSPWKL GDCGPQDDDTESCLCLEMN VAQEFQLRRRQLGSQGSSWS KWSSPVCVPPENPPQPQVRF SVEQLGQDGRRLTLKEQP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------