

Human IL-3RA/CD123 protein (Recombinant) (C-His) (STJP019904)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	3563
Gene Symbol	IL3RA
UniProt ID	IL3RA_HUMAN
Immunogen Region	Thr19-Arg305
Immunogen Sequence	TKEDPNPPITNLRMKAKAQQ LTWDLNRNVTDIECVKDADY SMPAVNNSYCQFGAISLCEV TNYTVRVANPPFSTWILFPE NSGKPWAGAENLTCWIHDVD FLSCSWAVGPGAPADVQYDL YLNVANRRQQYECLHYKTDA QGTRIGCRFDDISRLSSGSQ SSHILVRGRSAAFIPCTDK FVVFSQLIETPPNMTAKCN KTHSFMHWKMRSHFNRFKFRY ELQIQKRMQPVITEQVRDR
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.