

Human IL-2 R beta/CD122 protein (Recombinant) (C-His&Avi) (STJP019906)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	3560
Gene Symbol	IL2RB
UniProt ID	IL2RB_HUMAN
Immunogen Region	Ala27-Asp239
Immunogen Sequence	AVNGTSQFTCFYNSRANISC VWSQDGALQDTSCQVHAWPD RRRWNQTCELLPVSQASWAC NLILGAPDSQKLTTVDIVTL RVLCREGVRWRVMAIQDFKP FENLRMLAPISLQVVHVETH RCNISWEISQASHYFERHLE FEARTLSPGHTWEEAPLLTL KQKQEWICLETLPDTQYEF QVRVKPLQGEFTTWSPWSQP LAFRTKPAALGKD
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.