

Human IL-2 R beta/CD122 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019907)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3560
Gene Symbol	IL2RB
UniProt ID	IL2RB_HUMAN
Immunogen Region	Ala27-Asp239
Immunogen Sequence	AVNGTSQFTCFYNSRANISC VWSQDGALQDTSCQVHAWPD RRRWNQTCELLPVSQASWAC NLILGAPDSQKLTTVDIVTL RVLCREGVRWRVMAIQDFKP FENLRMLAPISLQVVHVETH RCNISWEISQASHYFERHLE FEARTLSPGHTWEEAPLLTL KQKQEWICLETLPDTQYEF QVRVKPLQGEFTTWSPWSQP LAFRTKPAALGKD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---