

Human IL-2 R alpha/CD25 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019910)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3559
Gene Symbol	IL2RA
UniProt ID	IL2RA_HUMAN
Immunogen Region	Glu22-Cys213
Immunogen Sequence	ELCDDDPPEIPHATFKAMAY KEGTMLNCECKRGFRRIKSG SLYMLCTGNSSHSSWDNQCQ CTSSATRNNTTKQVTPQPEEQ KERKTTEMQSPMQPVDQASL PGHCREPPPWENEATERIYH FVVGQMVYYQCVQGYRALHR GPAESVCKMTHGKTRWTQPQ LICTGEMETSQFPGEEKPQA SPEGRPESETSC

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---