

Human IL-1R1/CD121a protein (Recombinant) (No Tag) (STJP019911)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3554
Gene Symbol	IL1R1
UniProt ID	IL1R1_HUMAN
Immunogen Region	Leu18-Thr332
Immunogen Sequence	LEADKCKEREKIIIVSSAN EIDVRPCPLNPNEHKGITW YKDDSKTPVSTEQASRIHQH KEKLWFVPAKVEDSGHYCYV VRNSSYCLRIKISAKFVENE PNLCYNAQAIFKQKLPVAGD GGLVCPYMEFFKNENNELPK LQWYKDKPILLDNIHFSGV KDRLIVMNVAEKHRGNYTCH ASYTYLGKQYPITRVIEFIT LEENKPTRPVIVSPANETME VDLGSQLICNVTGQLSD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------