

Human IL-1 beta protein (Recombinant) (C-His&Avi) {Biotin} (STJP019913)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3553
Gene Symbol	IL1B
UniProt ID	IL1B_HUMAN
Immunogen Region	Ala117-Ser269
Immunogen Sequence	APVRSLNCTLRDSQQKSLVM SGPYELKALHLQGQDMEQQV VFMSMSFVQGEESNDKIPVAL GLKEKNLYLSCVLKDDKPTL QLESVDPKNYPKKKMEKRFV FNKIEINNKLFEFESAQFPNW YISTSQAENMPVFLGGTKGG QDITDFTMQFVSS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---