

Human FLT-3/FLK-2/CD135 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019916)

GENERAL INFORMATION

| | |
|---------------|-----------|
| Product Type | Proteins |
| Host / Source | CHO cells |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Supplied as 0.22 µg/ml filtered solution in PBS (pH 7.4). |
| Storage Instruction | Store at -70 °C. This product is stable at ≤-70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|---|
| Gene ID | 2322 |
| Gene Symbol | FLT3 |
| UniProt ID | FLT3_HUMAN |
| Immunogen Region | Asn27-Asn541 |
| Immunogen Sequence | NQDLPIKCVLINHKNDSS VGKSSSYPMVSESPEDLGCA LRPQSSGTVYEAHAVEVDVS ASITLQVLVDAPGNISCLWV FKHSSLNCQPHFDLQNRGVV SMVILKMTETQAGEYLLFIQ SEATNYTILFTVSIRNTLLY TLRRPYFRKMENQDALVCIS ESVPEPIVEWVLCDSQGESC KEESPAVVKKEEKVLHELFG TDIRCCARNELGRECTRLFT IDLNQTPQTTLPQLFLKVG |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|