

Human FGFR-3 beta (IIlb) /CD333 protein (Recombinant) (C-His&Avi) (STJP019918)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2261
Gene Symbol	FGFR3
UniProt ID	FGFR3_HUMAN
Immunogen Region	Asp127-Gly377
Immunogen Sequence	DAPSSGDDEDEGEDEAEDTGV DTGAPYWTRPERMDKKLLAV PAANTVRFRCPAAGNPTPSI SWLKNGREFRGEHRIGGIKL RHQQWSLVMESVVPDRGNY TCVVENKFGSIRQTYTLDVL ERSPHRPILQAGLPANQTAV LGSDVEFHCKVYSDAQPHIQ WLKHVEVNGSKVGPDPGTPYV TVLKSWISESVEADVRLRLA NVSERDGGEYLCRATNFIGV AEKAFWLSVHGPRAAEEEL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---