

## Human FGFR-3 alpha (IIlc) /CD333 protein (Recombinant) (C-His&Avi) (STJP019924)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |                                                                                                                                   |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                    |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                     |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">2261</a>                                                                                                                                                                                                                                                |
| Gene Symbol        | <a href="#">FGFR3</a>                                                                                                                                                                                                                                               |
| UniProt ID         | <a href="#">FGFR3_HUMAN</a>                                                                                                                                                                                                                                         |
| Immunogen Region   | Glu23-Gly375                                                                                                                                                                                                                                                        |
| Immunogen Sequence | ESLGTEQRVVGRAAEVP GPE PGQEQQLVFGSGDAVELSCP PPGGGPMGPTVWVKDGTGLV<br>PSERVLVGPQRLQVLNASHE DSGAYSCRQRLTQRVLCHFS VRVTDAPSSGDDEDEGEDEAE<br>DTGVDTGAPYWTRPERMDKK LLAVPAANTVRFRCPAAGNP TPSISWLKNGREFRGEHRIG<br>GIKLRHQQWLSLVMEVPSD RGNVTCVVENKFGSIRQTYT LDVLESPHRPILQAGLPA |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|