

Human FGFR-3 alpha (IIlc) /CD333 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019925)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2261
Gene Symbol	FGFR3
UniProt ID	FGFR3_HUMAN
Immunogen Region	Glu23-Gly375
Immunogen Sequence	ESLGTEQRVVGRAAEVP GPE PGQEQQLVFGSGDAVELSCP PPGGGPMGPTVWVKDGTGLV PSERVLVGPQRLQVLNASHE DSGAYSCRQRLTQRVLCHFS VRVTDAPSSGDDEDEGEDEAE DTGVDTGAPYWTRPERMDKK LLAVPAANTVRFRCPAAGNP TPSISWLKNGREFRGEHRIG GIKLRHQQWLSLVMEVPSD RGNVTCVVENKFGSIRQTYT LDVLESPHPRILQAGLPA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---