

## Human B7-H3 (4Ig) /B7-H3b/CD276 protein (Recombinant) (C-hFc) (STJP019927)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">80381</a>
Gene Symbol	<a href="#">CD276</a>
UniProt ID	<a href="#">CD276_HUMAN</a>
Immunogen Region	Gly27-Thr461
Immunogen Sequence	GALEVQVPEDPVVALVGTDA TLCCSFSPPEPGFSLAQLNLI WQLTDTKQLVHSFAEGQDQG SAYANRTALFPDLLAQGNAS LRLQRVRVADEGSFTCFVSI RDFGSAAVSLQVAAAPYSKPS MTLEPNKDLRPGDVTITCS SYQGYPEAEVFWQDGGQGVPL TGNVTTSQMANEQGLFDVHS ILRVVLGANGTYSCLVRNPV LQQDAHSSVTITPQRSPTGA VEVQVPEDPVVALVGTDAT

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---