

Human FcRn/FCGRT&B2M protein (Recombinant) (C-His&Avi) {Biotin} (STJP019929)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2217
Gene Symbol	FCGRT
UniProt ID	FCGRN_HUMAN
Immunogen Region	Ala24-Ser297
Immunogen Sequence	AESHLISLLYHLTAVSSPAPG TPAFWVSGWLGPQQYLSYNS LRGEAEPGAWVWENQVSWY WEKETDDLRIKEKLFLEAFK ALGGKGPYTLQGGLGCELGP DNTSVPTAKFALNGEEFMNF DLKQGTWGGDWPEALAIQR WQQQDKAANKELTFLLFSCP HRLREHLERGRGNLEWKEPP SMRLKARPSSPGFSVLTCSA FSFYPELQLRFLRNGLAAG TGQGDGFPNSDGSFHASSS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---