

Human Fc-gamma RIII beta/FCGR3B/CD16b (S36R, N82D, I106V) protein (Recombinant) (C-His) (STJP019933)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	2215
Gene Symbol	FCGR3B
UniProt ID	FCG3B_HUMAN
Immunogen Region	Gly17-Ser200
Immunogen Sequence	GMRTEDLPKAVVFLEPQWYS VLEKDSVTLKCGAYSPEDN STQWFHNENLISSQASSYFI DAATVNDSGEYRCQTNLSTL SDPVQLEVHIGWLLQAPRW VFKEEDPIHLRCHSWKNTAL HKVTYLQNGKDRKYFHHNSD FHIPKATLKDSGSYFCRGLV GSKNVSSETVNITITQGLAV STIS
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.