

Human Fc-gamma RIII alpha/FCGR3A/CD16a (F176V) protein (Recombinant) (C-His&Avi) {Biotin} (STJP019936)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	2214
Gene Symbol	FCGR3A
UniProt ID	FCG3A_HUMAN
Immunogen Region	Gly17-Gln208
Immunogen Sequence	GMRTEDLPKAVVFLEPQWYR VLEKDSVTLKCQGAYSPEDN STQWFHNESLISSQASSYFI DAATVDDSGEYRCQTNLSTL SDPVQLEVHIGWLLQAPRW VFKEEDPIHLRCHSWKNTAL HKVTYLQNGKGRKYFHHNSD FYIPKATLKDSGSYFCRGLV GSKNVSSETVNITITQGLAV STISSFFPPGYQ
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.