

Human Fc gamma RIIB/FCGR2B/CD32b protein (Recombinant) (C-His&Avi) (STJP019937)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2213
Gene Symbol	FCGR2B
UniProt ID	FCG2B_HUMAN
Immunogen Region	Ala46-Pro217
Immunogen Sequence	APPKAVLKLEPQWINVLQED SVTLTCRGTHSPESDSIQWF HNGNLIPHTHTQPSYRFKANN NDSGEYTCQTGQTSLSDPVH LTVLSEWLVLQTPHLEFQEG ETIVLRCHSWKDKPLVKVTF FQNGKSKKFSRSDPNFSIPQ ANHSHSGDYHCTGNIGYTLY SSKPVTITVQAP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---