

Human Fc gamma RIIA/FCGR2A/CD32a (H167R) protein (Recombinant) (C-His&Avi) {Biotin} (STJP019940)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	2212
Gene Symbol	FCGR2A
UniProt ID	FCG2A_HUMAN
Immunogen Region	Ala36-Ile218
Immunogen Sequence	AAPPKAVLKLEPPWINVLQE DSVTLTCQGARSPESDSIQW FHNGNLIPTHTQPSYRFKAN NNDSGEYTCQTGQTSLSDPV HLTVLSEWLVLQTPHLEFQE GETIMLRCHSWKDKPLVKVT FFQNGKSQKFSRLDPTFSIP QANHSHSGDYHCTGNIGYTL FSSKPVTTITVQVPSMGSSSP MGI
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.