

Human ErbB-2/HER2/CD340 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019947)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2064
Gene Symbol	ERBB2
UniProt ID	ERBB2_HUMAN
Immunogen Region	Thr23-Thr652
Immunogen Sequence	TQVCTGTDMLRLPASPETH LDMLRHLYQGCQVVQGNLEL TYLPTNASLSFLQDIQEVQG YVLIAHNQVRQVPLQRLRIV RGTQLFEDNYALAVLDNGDP LNNTTPVTGASPGGLRELQL RSLTEILKGGVLIQRNPQLC YQDTILWKDIFHKNNQLALT LIDTNRSRACHPCSPMCKGS RCWGESSEDCQSLTRTV CAG GCARCKGPLPTDCCHEQCAA GCTGPKHSDCLACLHFNHS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---