

## Human ErbB-2/HER2/CD340 protein (Recombinant) (C-hFc) (STJP019948)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2064</a>
Gene Symbol	<a href="#">ERBB2</a>
UniProt ID	<a href="#">ERBB2_HUMAN</a>
Immunogen Region	Thr23-Thr652
Immunogen Sequence	TQVCTGTDMLRRLPASPETH LDMLRHLYQGCQVVQGNLEL TYLPTNASLSFLQDIQEVQG YVLIAHNQVRQVPLQRLRIV RGTQLFEDNYALAVLDNGDP LNNTTPVTGASPGGLRELQL RSLTEILKGGVLIQRNPQLC YQDTILWKDIFHKNNQLALT LIDTNRSRACHPCSPMCKGS RCWGESSEDCQSLTRTV CAG GCARCKGPLPTDCCHEQCAA GCTGPKHSDCLACLHFNHS

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---