

Human HER1/ERBB1/EGFR protein (FITC-Labeled Recombinant) (C-His&Avi) (STJP019951)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1956
Gene Symbol	EGFR
UniProt ID	EGFR_HUMAN
Immunogen Region	Leu25-Ser645
Immunogen Sequence	LEEKKVCQGTSNKLTQLGTF EDHFLSLQRMFNNECVVLGN LEITYVQRNYDLSFLKTIQE VAGYVLIALNTVERIPLNL QIIRGNMYYENSALAVLSN YDANKTGLKELPMRNLQEIL HGAVRFSNNPALCNVESIQW RDIVSSDFLSNMSMDFQNLH GSCQKCDPSCPNGSCWGAGE ENCQKLTIIICAQQCSGRGR GKSPSDCCHNQCAAGCTGPR ESDCLVCRKFRDEATCKDT
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.