

Human HER1/ERBB1/EGFR (25-378) protein (Recombinant) (C-His&Avi) (STJP019952)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1956
Gene Symbol	EGFR
UniProt ID	EGFR_HUMAN
Immunogen Region	Leu25-Ser378
Immunogen Sequence	LEEKKGNYVVTDHGSCVRAC GADSYEMEEDGVRKCKKCEG PCRKVCNGIGIGEFKDSLSI NATNIKHFKNCTSIISGDLHI LPVAFRGDSFTHTPPLDPQE LDILKTVKEITGFLLIQAWP ENRTDLHAFENLEIIRGRK QHGQFSLAVVSLNITSLGLR SLKEISDGDVVISGNKNLCY ANTINWKKLFGTSGQKTKII SNRGENSCKATGQVCHALCS PEGCWGPEPRDCVSCRNV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---