

Human HER1/ERBB1/EGFR (25-378) protein (PE-Labeled Recombinant) (C-His&Avi) (STJP019955)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as 0.22 Mu m filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -70 °C. This product is stable at ≤ -70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1956
Gene Symbol	EGFR
UniProt ID	EGFR_HUMAN
Immunogen Region	Leu25-Ser378
Immunogen Sequence	LEEKKGNYVVTDHGSCVRAC GADSYEMEEDGVRKCKKCEG PCRKVCNGIGIGEFKDSL NATNIKHFKNCTSIISGDLHI LPVAFRGDSFTHTPLDPQE LDILKTVKEITGFLLIQAWP ENRTDLHAFENLEIIRGRK QHGQFSLAVVSLNITSLGLR SLKEISDGDVVIISGNKNCY ANTINWKKLFGTSGQKTKII SNRGENSCKATGQVCHALCS PEGCWWGPEPRDCVSCRNV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---