

Human CTLA-4/CD152 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019957)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1493
Gene Symbol	CTLA4
UniProt ID	CTLA4_HUMAN
Immunogen Region	Lys36-Asp161
Immunogen Sequence	KAMHVAQPAVVLASSRGIAS FVCEYASPGKATEVRVTVL R QADSQVTEVCAATYMMGNEL TFLDDSICTGTSSGNQVNLT IQGLRAMDTGLYICKVELMY PPPYYLGIGNGTQIYVIDPE PCPDSD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---