

## Human CSF2RA/GM-CSF R alpha/CD116 protein (Recombinant) (C-His) (STJP019958)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1438</a>
Gene Symbol	<a href="#">CSF2RA</a>
UniProt ID	<a href="#">CSF2R_HUMAN</a>
Immunogen Region	Glu23-Gly320
Immunogen Sequence	EKSDLRTVAPASSLNVRFDSTRMNLSDWDCQENTTFQSKCFLTDKKNRVVEPRLSNNECSCT FREICLHEGVTFEVHVNTSQ RGFQQKLLYPNSGREGTAAQ NFSCFIYNADLMNCTWARGP TAPRDVQYFLYIRNSKRRE IRCPYYIQDSGTHVGCHLDN LSGLTSRNYFLVNGTSREIG IQFFDSLDDTKIERFNPPS NVTVRCNTTHCLVRWKQPRT YQKLSYLDQYQLDVHRKN

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---