

## Human CSF2RA/GM-CSF R alpha/CD116 protein (Recombinant) (C-hFc) (STJP019959)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Formulation         | Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                    |

### TARGET INFORMATION

|                    |   |
|--------------------|---|
| Gene ID            | <a href="#">1438</a>  |
| Gene Symbol        | <a href="#">CSF2RA</a>  |
| UniProt ID         | <a href="#">CSF2R_HUMAN</a>   |
| Immunogen Region   | Glu23-Gly320  |
| Immunogen Sequence | EKSDLRTVAPASSLNVRFDSTRMNLSDWDCQENTTFFSKCFLETDKKNRVVEPRLSNNECSCT<br>FREICLHEGVTFEVHVNTSQ RGFQQKLLYPNSGREGTAAQ NFSCFIYNADLMNCTWARGP<br>TAPRDVQYFLYIRNSKRRE IRCPYYIQDSGTHVGCHLDN LSGLTSRNYFLVNGTSREIG<br>IQFFDSLDDTKIERFNPPS NVTVRCNTTHCLVRWKQPRT YQKLSYLDQYQLDVHRKN |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|