

Human CSF2RA/GM-CSF R alpha/CD116 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019960)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1438
Gene Symbol	CSF2RA
UniProt ID	CSF2R_HUMAN
Immunogen Region	Glu23-Gly320
Immunogen Sequence	EKSDLRTVAPASSLNVRFDSTMTNLSWDCQENTTFSKCFLTDKKNRVVEPRLSNNECSCT FREICLHEGVTFEVHVNTSQRGFQQKLLYPNSGREGTAAQNFSCFIYNADLMNCTWARGP TAPRDVQYFLYIRNSKRREIRCPYYIQDSGTHVGCHLDNLSGLTSRNYFLVNGTSREIG IQFFDSLDDTKIERFNPPSNVTVRCNTTHCLVRWKQPRTYQKLSYLDYQYQLDVRHN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---