

Human CEACAM8/CD66b protein (Recombinant) (C-His&Avi) {Biotin} (STJP019961)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1088
Gene Symbol	CEACAM8
UniProt ID	CEAM8_HUMAN
Immunogen Region	Gln35-Ser319
Immunogen Sequence	QLTIEAVPSNAAEGKEVLLL VHNLPQDPRGYNWYKGETVD ANRRIIGYVISNQQITPGPA YSNRETIYPNASLLMRNVTR NDTGSYTLQVIKLNLMSEEV TGQFSVHPETPKPSSISNNS NPVEDKDAVAFTCEPETQNT TYLWWWNGQSLPVSPRLQLS NGNRRTLTLVTRNDVGPYE CEIQNPASANFSDPVTNLVL YGPDAPTISPSTYYHAGVN LNLSCHAASNPPSQYSWSV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---