

Human CEACAM5/CEA/CD66e (145-322) protein (Recombinant) (C-His) (STJP019964)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as 0.22 µg/ml filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -70 °C. This product is stable at ≤-70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1048
Gene Symbol	CEACAM5
UniProt ID	CEAM5_HUMAN
Immunogen Region	Pro145-Pro322
Immunogen Sequence	PKPSISSNNSKPVEDKDAVA FTCEPETQDATYLWWVNNQS LPVSPRLQLSNGNRTLTLFN VTRNDTASYKCETQNPVSAR RSDSVILNVLYGPDAPTISP LNTSYRSGENLNLSCAASN PPAQYSWFVNGTFQQSTQEL FIPNITVNNNSGSYTCQAHNS DTGLNRTTVTITVYAEF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---