

Human CEACAM5/CEA/CD66e protein (Recombinant) (C-His) (STJP019966)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1048
Gene Symbol	CEACAM5
UniProt ID	CEAM5_HUMAN
Immunogen Region	Lys35-Ala685
Immunogen Sequence	KLTIESTPFNVAEGKEVLLL VHNLPQHLEFGYSWYKGERVD GNRQIIGYVIGTQQATPGPA YSGREIIPNASLLIQNIIQ NDTGFYTLHVIKSDLVNEEA TGQFRVYPELPKPSISSNNS KPVEDKDAVAFTCEPETQDA TYLWWWNNQSLPVSPRLQLS NGNRTLTLFNVTRNDTASYK CETQNPVSARRSDSVILNVL YGPDAPTISPLNTSYRSGEN LNLSCHAASNPPAQYSWFV
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.