

Human CEACAM5/CEA/CD66e protein (Recombinant) (C-His&Avi) {Biotin} (STJP019967)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 M µm filtered solution in 20 mM PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1048
Gene Symbol	CEACAM5
UniProt ID	CEAM5_HUMAN
Immunogen Region	Lys35-Ala685
Immunogen Sequence	KLTISTPFNVAEGKEVLLL VHNLPQHLEFGYSWYKGERVD GNRQIIGYVIGTQQATPGPA YSGREIIPNASLLIQNIQ NDTGFYTLHVIKSDLVNEEA TGQFRVYPELPKPSISSNNS KPVEDKDAVAFTCEPETQDA TYLWWWNNQSLPVSPRLQLS NGNRTLTLFNVTRNDTASYK CETQNPVSARRSDSVILNVL YGPDAPTISPLNTSYRSGEN LNLSCHAASNPPAQYSWFV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---