

## Human CD44 protein (Recombinant) (C-hFc) (STJP019977)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">960</a>
Gene Symbol	<a href="#">CD44</a>
UniProt ID	<a href="#">CD44_HUMAN</a>
Immunogen Region	Gln21-Pro220 & Gln386-Ala427
Immunogen Sequence	QIDLNITCRFAGVFHVEKNG RYSISRTEAADLCKAFNSTL PTMAQMEKALSIGFETCRYG FIEGHVVIPRIHPNSICAA N TG VYILTSNTSQYD TYCFN ASAPPEEDCTSVTDL PNAFD GPITIVNRD GTRYVQKGE YRTNPEDIYPSNPTDDDVSS GSSSERSSTSGGYFYTFST VHIPDEDESPWITDSTDRIP QATPSSTTEETATQKEQWFG NRWHEGYRQTPKEDSHSTT

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---