

Human CD44 protein (Recombinant) (C-His) (STJP019978)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
UniProt ID	CD44_HUMAN
Immunogen Region	Gln21-Pro220 & Gln386-Ala427
Immunogen Sequence	QIDLNITCRFAGVFHVEKNG RYSISRTEAADLCKAFNSTL PTMAQMEKALSIGFETCRYG FIEGHVVIPRIHPNSICAA NNTGVYILTSNTSQYDTCFN ASAPPEEDCTSVTDLPNAFD GPITITIVNRDGTTRYVQKGE YRTNPEDIYPSNPTDDDVSS GSSSERSSTSGGYIFYTFST VHPIPEDDSPWITDSTDRIP QATPSSTTEETATQKEQWFG NRWHEGYRQTPKEDSHSTT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---