

Human B7-H3 (4Ig) /B7-H3b/CD276 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019984)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80381
Gene Symbol	CD276
UniProt ID	CD276_HUMAN
Immunogen Region	Gly27-Thr461
Immunogen Sequence	GALEVQVPEDPVVALVGTDA TLCCSFSPPEPGFSLAQLNLI WQLTDTKQLVHSFAEGQDQG SAYANRTALFPDLLAQGNAS LRLQRVRVADEGSFTCFVSI RDFGSAAVSLQVAAAPYSKPS MTLEPNKDLRPGDVTITCS SYQGYPEAEVFWQDGGQGVPL TGNVTTSQMANEQGLFDVHS ILRVVLGANGTYSCLVRNPV LQQDAHSSVTITPQRSPTGA VEVQVPEDPVVALVGTDAT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---