

Human TNFRSF5/CD40 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019986)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	958
Gene Symbol	CD40
UniProt ID	TNRS_HUMAN
Immunogen Region	Glu21-Arg193
Immunogen Sequence	EPPTACREKQYLINSQCCSL CQPGQKLVSDCTEFTETECL PCGESEFLDTWNRETHCHQH KYCDPNLGLRVQKGTSETD TICTCEEGWHCTSEACESCV LHRSCSPGFGVKQIATGVSD TICEPCPVGFFSNVSSAFEK CHPWTSCETKDLVVQQAGTN KTDVVCGPQDRLR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---