

Human CD38 protein (Recombinant) (C-His&Avi) {Biotin} (STJP019988)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	952
Gene Symbol	CD38
UniProt ID	CD38_HUMAN
Immunogen Region	Val43-Ile300
Immunogen Sequence	VPRWRQQWSGPGTTKRFPEI VLARCVKYTEIHPEMRHVDC QSVWDAFKGAFISKHPCNIT EEDYQPLMKLGTQTVPCNKI LLWSRIKDLAQFTQVQRDM FTLEDTLGLYADDLTWCGE FNTSKINYQSCPDWRKDCSN NPVSVFWKTVSRRFAEAACD VVHVMLNGSRSKIFDKNSTF GSVEVHNLQPEKVQTLAEWV IHGGREDSRDLCDPTIKEL ESIISKRNQFSCKNIRP
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.