

Human Siglec-3/CD33 protein (Recombinant) (C-His&Avi) (STJP019990)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in 20 mM MES, 500 mM NaCl, 200 mM L-Arginine (pH 5.0). Normally 8% mannitol is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	945
Gene Symbol	CD33
UniProt ID	CD33_HUMAN
Immunogen Region	Asp18-His259
Immunogen Sequence	DPNFWLQVQESVTVQEGLCV LVPCTFFHPIPYDDKNSPVH GYWFREGAIISRDSPVATNK LDQEVQEETQGRFRLLGDPS RNNCSLSIVDARRRDNGSYF FRMERGSTKYSPQLSVH VTDLTHRPKILIPGTLEPGH SKNLTCSVSWACEQGTPIF SWLSAAPTSLGPRTHSSVL IITPRPDHGTNLTCQVKFA GAGVTTERTIQLNVTYVPQN PTTGIFPGDGSGKQETRAG
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.