

Human TNFRSF8/CD30 protein (Recombinant) (C-hFc) (STJP019995)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	943
Gene Symbol	TNFRSF8
UniProt ID	TNR8_HUMAN
Immunogen Region	Phe19-Lys379
Immunogen Sequence	FPQDRPFEDTCHGNPSHYD KAVRRCCYRCPMGLFPTQQC PQRPTDCRKQCEPDYYLDEA DRCTACVTCSRDDLVEKTPC AWNSSRVCECRPGMFCSTSA VNSCARCFHSHVCPAGMIVK FPGTAQKNTVCEPASPGVSP ACASPENCKEPSSGTIPQAK PTPVSPATSSASTMPVRGGT RLAQEAASKLTRAPDSPSSV GRPSSDPGLSPTQPCPEGSG DCRKQCEPDYYLDEAGRCT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---