

Human B7-2/CD86 protein (Recombinant) (C-hFc) (STJP019996)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	942
Gene Symbol	CD86
UniProt ID	CD86_HUMAN
Immunogen Region	Leu26-Pro247
Immunogen Sequence	LKIQAYFNETADLPCQFANS QNQSLSELVVFWDQENLVL NEVYLGKEKFDSVHSKYMGR TSFSDSDSWTLRLHNLQIKDK GLYQCIIHHKKPTGMIRIHQ MNSELSVLANFSQPEIVPIS NITENVYINLTCSSIHGYPE PKKMSVLLRTKNSTIEYDGV MQKSQDNVTELYDVSISLSV SFPDVTSNMTIFCILETDKT RLLSSPFSIELEDPQPPPDH IP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---