

Human/Cynomolgus/Rhesus macaque CD28 protein (Recombinant) (N-mFc) (STJP020002)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	940
Gene Symbol	CD28
UniProt ID	CD28_HUMAN
Immunogen Region	Asn19-Pro152
Immunogen Sequence	NKILVKQSPMLVAYDNAVNLSCKYSYNLFSREFRASLHKG LDSAVEVCVVYGNYSQQLQV YSKTGFNCDBGKLGNESVTFY LQNLVYNQTDIYFCKIEVMY PPPYLDNEKSNGTIIHVKGK HLCPSPLFPGPSKP
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.