

Human CD3E&CD3G protein (Recombinant) (C-hFc) (STJP020009)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Asp23-Asp126
Immunogen Sequence	DGNEEMGGITQTPYKVSISG TTVILTCPQYPGSEILWQHN DKNIGGDEDDKNIGSDEDHL SLKEFSELEQSGYYVCYPRG SKPEDANFYLYLRARVCENC MEMDQSIKGNHLVKVYDYQE DGSVLLTCDAEAKNITWFKD GKMIGFLTEDKKKWNLGSNA KDPRGMYQCKGSQNKSKPLQ VYYRMCQCNCIELNAATIS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---