

Human MS4A1/CD20 protein (Recombinant) (C-His&Avi) (STJP020011)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as 0.22 µg/ml filtered solution in 20 mM PB, 10% Glycerol, 0.4 M Urea, 1 mM GSH, 0.1 mM GSSG (pH 7.0).
Storage Instruction	Store at -70 °C. This product is stable at ≤ -70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	931
Gene Symbol	MS4A1
UniProt ID	CD20_HUMAN
Immunogen Region	Ile141-Ser188
Immunogen Sequence	IKISHFLKMESLNFIRAHTP YINIYNCEPANPSEKNSPST QYCYSIQS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---