

Human CD7 protein (Recombinant) (C-His&Avi) {Biotin} (STJP020014)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	924
Gene Symbol	CD7
UniProt ID	CD7_HUMAN
Immunogen Region	Ala26-Pro180
Immunogen Sequence	AQEVQQSPHCTTVPGASVN ITCSTSGGLRGIYLRQLGPQ PQDIIYYEDGVVPTTDRRFR GRIDFSGSQDNLITMHRLO LSDTGTYTCQAITEVNVYGS GTLVLVTEEQSQGWHRCSDA PPRASALPAPPTGSALPDPQ TASALPDPPAASALP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---