

Human CD7 protein (PE-Labeled Recombinant) (C-His&Avi) (STJP020015)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as 0.22 µg/ml filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -70 °C. This product is stable at ≤ -70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	924
Gene Symbol	CD7
UniProt ID	CD7_HUMAN
Immunogen Region	Ala26-Pro180
Immunogen Sequence	AQEVQQSPHCTTVPVGASVN ITCSTSGGLRGIYLRQLGPQ PQDIIYYEDGVVPTDRRFR GRIDFSGSQDNLITMHRQL LSDTGTYTCAITEVNVYGS GTLVLVTEEQSQGWHRCSDA PPRASALPAPPTGSALPDPQ TASALPDPPAASALP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---