

Human CD5 protein (Recombinant) (C-hFc&Avi) {Biotin} (STJP020019)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	921
Gene Symbol	CD5
UniProt ID	CD5_HUMAN
Immunogen Region	Arg25-Asn371
Immunogen Sequence	RLSWYDPDFQARLTRNSKQC QGQLEVLKDGWHMVCSQSW GRSSKQWEDPSQASKVCQRL NCGVPLSLGPFLLVITYTPQSS IICYGQLGSFNSCHSRNDM CHSLGLTCLPQKTTPTTR PPPTTTPEPTAPPRLQLVAQ SGGQHCAGVVEFYSGSLGGT ISYEAQDKTQDLENFLC>NNL QCGSFLKHLPETEAGRAQDP GEPREHQPLPIQWKIQNSSC TSLEHCFRRIKIPQKSGRVL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---