

Human LEU3/CD4 protein (Recombinant) (C-His&Avi) {Biotin} (STJP020020)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	920
Gene Symbol	CD4
UniProt ID	CD4_HUMAN
Immunogen Region	Lys26-Trp390
Immunogen Sequence	KKVVLGKKGDTVELTCTASQ KKSIIQFHWKNSNIKILGNQ GSFLTGKPSKLNDRADSRRS LWDQGNFPLIINKLKIEDSD TYICEVEDQKEEVQLLVFGL TANSDTHLLQGQSLTLTLES PPGSSPSVQCRSPRGKNIQG GKTLVSVSQLELQDSGTWTCT VLQNQKKVEFKIDIVVLAFAQ KASSIVYKKEGEQVEFSFPL AFTVEKLTGSGELWWQAERA SSSKSWITFDLKNKEVSVK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---