

Human TNFRSF17/BCMA/CD269 Trimer protein (Recombinant) (C-His&Avi) {Biotin} (STJP020030)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	608
Gene Symbol	TNFRSF17
UniProt ID	TNR17_HUMAN
Immunogen Region	Met1-Ala54
Immunogen Sequence	MLQMAGQCSQNEYFDSLLHA CIPCQLRCSSNTPPLTCQRY CNASVTNSVKGTNA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---