

## Human ALCAM/CD166 protein (Recombinant) (C-His&Avi) {Biotin} (STJP020034)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">214</a>
Gene Symbol	<a href="#">ALCAM</a>
UniProt ID	<a href="#">CD166_HUMAN</a>
Immunogen Region	Trp28-Ala526
Immunogen Sequence	WYTVNSAYGDTIIIPCRLDV PQNLMFGKWKYEKPDGSPVF IAFRSSTKKSQYDDVPEYK DRLNLSENYTLISNARISD EKRFVCMVLVTEDNVFEPTI VKVFKQPSKPEIVSKALFLE TEQLKKGDCISEDSPDGN ITWYRNGKVLHPLEGAVVII FKKEMDPVTQLYTMTSTLEY KTTKADIQMPFTCSVTYYGP SGQKTIHSEQAVFDIYYPTE QVTIQVLPPKNAIKEGDNI

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------